



Boss Cryo-Hose™



When an injury occurs call on the Boss™. The new innovative Boss™ drinking system will provide the hydration your athletes needs as well as cryo-therapy when an injury occurs.

Using Aircast's® famous cryo-cuffs and the Boss™ cryo-attachment accessory, you can quickly provide cold



treatment on the sideline when needed. The treatment is simple and starts immediately. First, attach the cryo-hose to the Boss™ valve and quick-connect the selected body part Aircast® cuff to the system. Second, fill the Aircast® cuff to the desired pressure and release the hose. When treatment is finished connect

the cryo-hose to the Aircast® cuff and empty the water into a sink, shower, drain or on the ground. When time is critical in treating inflammation, look professional and provide the treatment that your athlete's deserve. The possibilities are endless when you have a Boss™ drinking system on your team.



Each Cryo/Cuff is anatomically designed to provide maximum coverage to the injured area. The **Knee, Shoulder, Ankle, Foot, and Thigh** Cryo/Cuffs are available in multiple sizes for expanded coverage. The Pediatric Ankle Cryo/Cuff and Pediatric Knee/Elbow Cryo/Cuff are specifically designed for children's injuries. Each Cryo/Cuff can be ordered as a single cuff or with the Boss™ Cryo-Attachment.

1-888-463-5699
www.outdoorboss.com
Copyright Go Flow, Inc.
2007